Document Date: June 10, 2016

Pasou Foods, Inc 1103 S. Huntington St. Syracuse, In 46567 Contact: Michele Pavel

Product: Hot Stix Item# 112436HS

Number of Servings 1 Serving Size : 1 Stix 55gm Case weight: 13 lbs

#### **Nutritional Facts Panel**

Nutrition Facts Serving Size Per container 100 Amount Per Serving 1 Calories 164	
Calories from Fat 70	Daily Value %
Total Fat 7.76 gm Saturated Fat 1.9 gm Trans Fat 0.3 gm *Poly unsaturated Fat 2.9 gm	12% 10%
*Monounsaturated Fat 2.7 gm Cholesterol 11 mg Sodium 370 mg Total Carbohydrate 17.66 gm Dietary Fiber 1.0 gm Sugars 1.8 gm Protein 6.00 gm Vitamin A < 10 IU Vitamin C < 0.1 IU Calcium 23 mg	3% 15% 6% 4%
Calcium 23 mg Iron 1.18 mg	

# Allergens

Milk Soy Wheat Yellow #5 & # 6 Wheat Gluten Sulfites

Wrapper may contain egg Keep Frozen or refrigerated Shelf Life 6 Months Product Description: Beef, Cheese, Seasoniongs in a Crispy wrapper.

## Ingredients

Water, Ground Beef, Textured Vegetable Protein(Soy flour, caramel coloring (contains Sulfites) Jalapeno Cheese (cheese whey, water, palm oil, modified food starch, cheddar cheese(milk cultures, salt enzymes,) modifiedtapioca starch, jalapeno peppers(water, aectic acid, salt, calcium chloride), maltodextrin, salt, sodium phosphate, sodium citrate, vinegar, sodium stearoyllactylate, natural flavor, mono and diglyceirides, sodium hexametaphosphate, annatto color, disodium inosinate, FD & C yellow 6. Jalapeno peppers(peppers, vinegar, salt, natural flavor, calcium chloride, sodium benzoate (preservative) sodium metabisulfite) Taco seasoning(dehydrated onion & garlic, salt spices, monosodium glutamate, red crushed pepper.wrapper ingredients Wheat flour, water, salt, wheat gluten, ascorbic acid,FD & C yellow #5 & #6 and cornstarch as an anti caking agent.

Case Weight: 13 lbs

# **Cooking Instructions**

## Deep Fryer

Recommend Thawing Preheat Oil 350 degrees Fry 3 to 4 minutes (Until Golden Brown)

### Bake:

Recommend Thawing Preheat oven to 450 degrees Bake 8 to 10 minutes

Approved by: Wayne E Glieber, Phd. Great lakes Scientific (616) 429-1000